- I - SISTEMA DI MISURA DELLA DIURESI ORARIA
E RELATIVI ACCESSORI: "MISURO"
- GB - SYSTEM FOR MEASURING HOURLY DIURESIS
AND RELEVANT ACCESSORIES: "MISURO"
- F - SYSTEME DE MESURE DE LA DIURESE HORAIRE
ET ACCESSOIRES ORRESPONDANTS: "MISURO"
- D - STUNDENDIURESE MESSYSTEM UND
ZUBEHÖR: "MISURO"
- E - SISTEMA DE MEDICIÓN DE LA DIURESIS
HORARIA Y ACCESORIOS CORRESPONDIENTES: MISURO
- P - SISTEMA PARA MEDIÇÃO DA DIURESE HORÁRIA E
RESPECTIVOS ACESSÓRIOS: "MISURO"
- N - ANORDNING TIL MÅLING AV EN TIMES DIURESE OG
TILBEHØR: "MISURO"
- NL - MEETSISTEEM VAN DE DIURESE PER UUR EN
BIJBEHORENDE ACCESSOIRES: "MISURO"
- FIN - TUNTIMÄÄRÄISEN VIRTSANERITYKSEN MITTAUSJÄRJESTELMÄ
JA SIIHEN KUULUVAT LISÄVARUSTEET: "MISURO"
- S - MÄTSYSTEM FÖR TIMDIURES OCH IFRÅGAVARANDE
TILLBEHÖR: "MISURO"
- DK - ANORDNING TIL MÅLING AF EN TIMES DIURESE
OG TILBEHØR: "MISURO"
- GR - ΣΥΣΤΗΜΑ ΜΕΤΡΗΣΗΣ ΤΗΣ ΩΡΙΑΙΑΣ ΔΙΟΥΡΗΣΗΣ
ΚΑΙ ΣΧΕΤΙΚΑ ΑΞΕΣΟΥΑΡ: "MISURO"
- PL - SYSTEM ZAMKNIĘTY DO POMIARU DIUREZY GODZINOWEJ

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APPLICATIONS

Measurement and collection of hourly diuresis.

INSTRUCTIONS FOR USE

- Thoroughly wash the hands and open the packet.
- Fix the device to the right-hand or left-hand side of the bed by means of the clamps (FIG 1A) or the hook (if present in the package) which must be attached to the handle of the device beforehand (FIG.1B). The bag must be fixed in a position not higher than the bladder of the patient.
- Put on a pair of disposable gloves (not provided).
- Remove the protective cap from the catheter cone (1) and connect it to the urological catheter (not provided) already inserted in the patient being careful not to contaminate it.
- Take off the gloves and discard them.
- Ensure that the tube (2) is positioned so as to allow regular urine flow. To do this use the clip (10) if present.
- When using bag models which have closing devices on the drainage tube (e.g. clamps, tube-clamp plates (3), roller clamps) bladder gymnastics can be carried out.

HOW TO DRAW A URINE SAMPLE

- Inform the patient and clamp the drainage tube beyond the drawing point.
- Thoroughly wash the hands.
- Put on a pair of disposable gloves (not provided).
- Disinfect the drawing point and insert the syringe needle (not provided) in the drawing point of the catheter cone.
- Aspirate the required quantity of urine and reopen the tube to restore drainage to the bag.

HOW TO MEASURE THE HOURLY DIURESIS AND EMPTY THE RIGID COLLECTION CONTAINER

- When the rigid container is in vertical position, the hourly diuresis can be measured using the graduated scale on the container (FIG 3).
- To empty the content of the rigid container into the collection bag, turn it anticlockwise (FIG 2) positioning it on the catch on the support handle in a horizontal position.
- Wait until the entire content has flowed into the bag.
- Reposition the rigid container in the vertical position (FIG 3) releasing it from the catch and turning it clockwise.

HOW TO REPLACE THE COLLECTION BAG.

CAUTION: The bag must be replaced after emptying out the rigid container.

- Thoroughly wash the hands.
- Put on a pair of disposable gloves (not provided).
- Position the container in the horizontal position (FIG 2).
- Close the tube-clamp plate (4) if present.
- Remove the drainage tube connected to the rigid container from its housing (5), initially exercising slight pressure in the direction of the container and then pulling in the opposite direction.
- Remove the drainage tube of the bag from the tube support of the handle (6).
- Remove the bag handle from its housing (7) exercising slight upward pressure.
- Open the packet containing the spare bag.
- Fix the handle of the spare bag in its housing (7) exercising slight pressure.
- Fix the tube of the spare bag inside the support (6) and connect it to the rigid measuring container in the position (5).
- Reposition the rigid container in the vertical position (FIG 3) releasing it from the catch and turning it clockwise.

HOW TO EMPTY THE BAG

Models with drainage tap :

- Thoroughly wash the hands and put on a pair of gloves (not provided).
 - Open the drainage tap and wait for the fluid in the bag to completely flow into the collection system adopted (container not provided).
 - When the operation has been completed, close the tap and disinfect it with gauze soaked in disinfectant (not provided) and replace it in its storage pocket (if present).
 - When the container has been emptied, wash and disinfect it. Take off the gloves.
- Models without drainage tap.
- For bags without tap, proceed with the bag replacement procedure described in the previous point.

HOW TO REPLACE "MISURO"

- Put on a pair of disposable gloves (not provided).
- Extremely gently disconnect the bag from the urological catheter of the patient, repositioning the protective cap on the catheter cone and dispose of the device.
- To apply a new "Misuro" on the patient, repeat the operations described previously in the INSTRUCTIONS FOR USE.

WARNINGS:

- For use only by qualified physicians. BIOSERVICE shall not be held responsible for patient injury in case of improper or incorrect use of the device, incorrect procedures, or in case of use by a non-qualified personnel.
- During handling of the device always use a sterile procedure in order to reduce the risks of infection.
- During use never fix the bag in a position higher than the bladder of the patient to prevent urine from flowing back into the bladder.
- Do not create traction on the catheter since this is a certain source of internal injury to the urethra.
- Do not use the bag for infusion.
- Do not fill the bag over its maximum capacity.
- The scale on the bag is expressed in millilitres (ml) and is approximate.
- Ensure that the tube is positioned so as to allow regular urine flow, thus preventing the urine from flowing back into the bladder.
- Do not place the bag on the floor when it is connected to the catheter, since this is a certain source of infection and malfunctioning of the bag.
- Only specially designed Bioservice spare bags may be connected to "Misuro".
- Maximum utilisation time: 15 days. Depending on the clinical picture of the patient, the medical/nursing staffs responsible for defining the need for more frequent replacement of the device.
- EO sterilized. Do not re-sterilize.
- Sterility is guaranteed if the packaging is undamaged.
- Disposable single-use/single-patient: re-use of the disposable may cause cross infection to patient and/or user.
- Pyrogen-free
- Expiry: 5 years with undamaged packaging and if stored in normal storage conditions (keep in a fresh and dry place, avoid the exposure to light and high temperature)
- To safeguard personnel, when required by hospital regulations or by the department, it might be necessary to wear sterile gloves, masks and protective glasses (Warning: these accessories are not supplied with the device).
- The urological catheter is not included in the product.
- Dispose of after use in conformity with the legislation in force.
- Do not use devices for premature infants, newborns and children under 10 kg.

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