- I - DEFLUSSORE PER LAVAGGIO UROLOGICO E ARTROSCOPICO
- GB - WASH SET FOR UROLOGICAL AND ARTHROSCOPIC WASHING
- F - CIRCUIT DE LAVAGE POUR LE RINCAGE UROLOGIQUE  
ET IRRIGATION ARTHROSCOPIQUE
- D - WASCHSYSTEM FÜR UROLOGISCHE UND ARTHROSKOPISCHE SPÜLUNG
- E - REGULADOR DE FLUJO PARA LAVADO UROLÓGICO Y ARTROSCÓPICO
- P - SET PARA LAVAGEM UROLÓGICA E ARTROSCÓPICA
- N - TILFØRSELSANORDNING FOR SKYLLING I FORBINDELSE MED UROLOGI  
OG ARTROSKOPI
- NL - SPOELINRICHTING VOOR UROLOGISCHE EN ARTROSCOPISCHE SPOELING
- FIN - VIRTAUSLAITE UROLOGISEEN JA ARTROSKOOPPISEEN PESUUN
- S - UROLOGISK OCH ARTROSKOPISK SKÖLJANORDNING
- DK - TILFØRSELSANORDNING TIL SKYLNING I FORBINDELSE MED UROLOGI  
OG ARTHROSKOPI

I

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GB

## APPLICATIONS

In urology it is used as connection system between the washing fluid container and the urological catheter or the cystoscope.

In arthroscopy it is used as connection system between the washing fluid container and the instrument for performing arthroscopy.

## INSTRUCTIONS FOR USE

- Thoroughly wash the hands and open the packet.
- Put on a pair of sterile gloves (not provided).
- Close the hose clamps (or the roller clamp for 1-way wash sets) on the line and connect the perforator(s) to the fluid container.
- Prime the tubes in such a way that no air remains in the circuit.
- Connect the end connector directly to the patient catheter or to the cystoscopy/arthroscopy instrument being careful not to contaminate it. The tube segment supplied in the package may be used to connect to the latter.
- Take off the gloves.
- Open the clamps and adjust the flow by means of the roller clamps.
- At the end of the procedure, wait until the tubes are empty and then close the clamps and roller clamps.
- Disconnect the end connector from the patient catheter or from the cystoscopy.
- Replace the protective caps of the perforators and the end connector.

## WARNINGS:

- Usage allowed only to qualified physicians. BIOSERVICE declines every responsibility for damages caused to the patient when the device is used in an improper way, with wrong manoeuvres or if handled by non qualified personnel.
- During handling of the device always use a sterile procedure in order to reduce the risks of infection.
- Do not create traction on the catheter since this is a certain source of internal injury to the urethra.
- Do not use to administer medication.
- Ensure that the tube is positioned so as to allow regular fluid flow.
- EtO sterilized. Do not re-sterilize.
- Sterility is guaranteed if the packaging is undamaged.
- Disposable single-use/single-patient: re-use of the disposable may cause cross infection to patient and/or user.
- Pyrogen-free
- Expiry: 5 years with undamaged packaging and if stored in normal storage conditions (keep in a fresh and dry place, avoid the exposure to light and high temperature).
- To safeguard personnel, when required by hospital regulations or by the department, it might be necessary to wear sterile gloves, masks and protective glasses (Warning: these accessories are not supplied with the device).
- The urological catheter is not included in the product.
- The washing solution is not included in the product.
- Dispose of after use in conformity with the legislation in force.
- Do not use devices for premature infants, newborns and children under 10 kg.

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