


CLAMPING DEVICES




- SMALL CLAMP ON/OFF FOR TUBES
- MEDIUM CLAMP ON/OFF FOR TUBES
- OPEN SLIDE CLAMP FOR TUBES
- CLAMP SHEET FOR BAG


SMALL CLAMP ON/OFF FOR TUBES

Code	Small clamp on/off	Material	Color	STD pack
C13203131 ETO - Radiation	 For tubes Ø 4.1 / Ø 5.5 mm	PP	White	8.000 pcs
C13203132 ETO - Radiation		PP	Blue	8.000 pcs
C13203133 ETO - Radiation		PP	Red	8.000 pcs


MEDIUM CLAMP ON/OFF FOR TUBES

Code	Medium clamp on/off	Material	Color	STD pack
C13503131 ETO - Radiation	 For tubes Ø 5.5 / Ø 6.8 mm	PP	White	3.000 pcs
C13503132 ETO - Radiation		PP	Blue	3.000 pcs
C13503133 ETO - Radiation		PP	Red	3.000 pcs

OPEN SLIDE CLAMP FOR TUBES

Code	Open slide clamp for tubes	Material	STD pack
B10003B ETO - Radiation		ABS	5.000 pcs

CLAMP SHEET FOR BAG

Code	Clamp sheet for bag	Material	STD pack
B10004B ETO - Radiation		ABS	4.000 pcs

● Dimension : «mm» ● Drawing scale: 1:1