



Components

- **FIXED MALE LUER LOCK**
- **FEMALE LUER LOCK**
- **MALE LUER LOCK CAPS**
- **PROTECTIVE CAP**
- **FIXED MALE LUER LOCK WITH VENTED CAP**
- **INJECTION SITES**

BioService Tunisie








Z.I Enfidha 4030 - Tunisia

Tel. +216 73 380 252 - Fax +216 73 380 288


e-mail : contact@bioservicetunisie.com

www.bioservicetunisie.com








FIXED MALE LUER LOCK

Code	MALE LUER LOCK	Material	Color	Quote	STD pack
C13304020 ETO - Radiation		ABS	Transparent	Ø 2.0 mm	20.000 pcs
C13314020 ETO - Radiation		ABS	Transparent	Ø 2.5 mm	20.000 pcs
C13374020 ETO - Radiation		ABS	Transparent	Ø 3.2 mm	20.000 pcs
C13384020 ETO - Radiation		ABS	Transparent	Ø 3.7 mm	20.000 pcs
C13324020 ETO - Radiation		ABS	Transparent	Ø 4.1 mm	20.000 pcs
C13344020 ETO - Radiation		ABS	Transparent	Ø 5.5 mm	20.000 pcs
C13354020 ETO - Radiation		ABS	Transparent	Ø 6.8 mm	20.000 pcs





PROTECTIVE CAP

Code	Vented cap for fixed male luer lock	Material	Color	Quote	STD pack
C16702130 ETO - Radiation (neutral)	 for item C133	LDPE	Neutral	Ø 9.6 mm	15.000 pcs








FEMALE LUER LOCK

Code	FEMALE LUER LOCK	Material	Color	Quote	STD pack
C11061030- ETO C11068010-Radiation		Rigid PVC	Transparent	Ø 2 mm	15.000 pcs
C11021030- ETO C11028010-Radiation		Rigid PVC	Transparent	Ø 2.5 mm	15.000 pcs
C11081030- ETO C11088010-Radiation		Rigid PVC	Transparent	Ø 3.2 mm	15.000 pcs
C11071030- ETO C11078010-Radiation		Rigid PVC	Transparent	Ø 3.7 mm	15.000 pcs
C11031030- ETO C11038010-Radiation		Rigid PVC	Transparent	Ø 4.1 mm	15.000 pcs
C11091030- ETO C11098010-Radiation		Rigid PVC	Transparent	Ø 5.0 mm	15.000 pcs
C11041030- ETO C11048010-Radiation		Rigid PVC	Transparent	Ø 5.5 mm	15.000 pcs



MALE LUER LOCK CAPS

Code	CLOSED CAP MALE LUER LOCK	Material	Color	STD pack
C10402120 ETO - Radiation		LDPE	Natural	30.000 pcs
C10402121 ETO - Radiation		LDPE	White	30.000 pcs
C10404121 ETO - Radiation		ABS	White	30.000 pcs
C10404145 ETO - Radiation		ABS	Green	30.000 pcs

FIXED MALE LUER LOCK WITH VENTED CAP

Code	MALE LUER LOCK WITH VENTED CAP	Material	Quote	STD pack
A01409990 ETO - Radiation		ABS-LDPE	Ø 2.0 mm	10.000 pcs
A01419990 ETO - Radiation		ABS-LDPE	Ø 2.5 mm	10.000 pcs
A01479990 ETO - Radiation		ABS-LDPE	Ø 3.2 mm	10.000 pcs
A01489990 ETO - Radiation		ABS-LDPE	Ø 3.7 mm	10.000 pcs
A01429990 ETO - Radiation		ABS-LDPE	Ø 4.1 mm	10.000 pcs
A01449990 ETO - Radiation		ABS-LDPE	Ø 5.5 mm	10.000 pcs
A01459990 ETO - Radiation		ABS-LDPE	Ø 6.8 mm	10.000 pcs

INJECTION SITES

Code	INJECTION SITE WITH MALE LUER LOCK	Material	Color	STD pack
A03009390 ETO - Radiation		ABS - Polysoprene (latex free)	Neutral	20.000 pcs
A01459990 ETO - Radiation		ABS - Polysoprene (latex free)	Yellow	20.000 pcs