


PVC TUBES FOR SURGICAL ASPIRATION



- Tube with wide lumen
- Light and easy to handle
- Non-deformable during aspiration
- Standard connectors



TUBES FOR SURGICAL ASPIRATION

The tubes are extremely soft and have a wide lumen to guarantee excellent flow of the aspirated fluid.

The tube has been designed in such a way as to prevent lumen constrictions during the aspiration phases.

Internally smooth and without unevenness to facilitate flow of the aspirated fluid. Externally the "satined" surface prevents undesired friction during use.

The suction tube is available in the versions:

- Female/Female (this version is fitted with a rigid Male/Male connector)
- Female/Male
- Female/Vacutip

1 ► The "Female" connector has a bellows geometry which allows easy handling and bending over 90° without reducing the lumen.

2 ► The "Male" connector, thanks to its rigidity, allows easy and secure fitting to any type of soft "Female" connector.

All the connectors guarantee perfect seal.

The device is supplied in double packaging sterilized by EtO.

On specific request of the customer the tubes can be supplied in non-sterile version.

SUCTION TUBES FEMALE/FEMALE CH25		
CODE	LENGTH	PACKAGING
H91510E	100 cm	50 pcs
H91515E	150 cm	50 pcs
H91518E	180 cm	50 pcs
H91520E	200 cm	50 pcs
H91522E	220 cm	50 pcs
H91525E	250 cm	50 pcs
H91527E	270 cm	40 pcs
H91530E	300 cm	40 pcs
H91535E	350 cm	30 pcs
H91537E	370 cm	30 pcs
H91540E	400 cm	30 pcs

SUCTION TUBES FEMALE/FEMALE CH30		
CODE	LENGTH	PACKAGING
H91380E	100 cm	50 pcs
H91385E	150 cm	50 pcs
H91388E	180 cm	50 pcs
H91390E	200 cm	50 pcs
H91392E	220 cm	50 pcs
H91395E	250 cm	40 pcs
H91397E	270 cm	40 pcs
H91400E	300 cm	40 pcs
H91405E	350 cm	30 pcs
H91407E	370 cm	30 pcs
H91410E	400 cm	30 pcs

SUCTION TUBES FEMALE/MALE CH25		
CODE	LENGTH	PACKAGING
H91270E	200 cm	50 pcs
H91280E	300 cm	40 pcs

SUCTION TUBES FEMALE/MALE CH30		
CODE	LENGTH	PACKAGING
H91330E	200 cm	50 pcs
H91340E	300 cm	40 pcs

SUCTION TUBES FEMALE/VACUTIP CH25		
CODE	LENGTH	PACKAGING
H91770E	200 cm	50 pcs
H91780E	300 cm	40 pcs

SUCTION TUBES FEMALE/VACUTIP CH30		
CODE	LENGTH	PACKAGING
H91810E	200 cm	50 pcs
H91820E	300 cm	40 pcs

NON-STERILE TUBING ON REEL

The tubing on reel gives the operators great versatility and economy in assembling customised circuits.

The dispenser packaging suitably protected the reel.

NON-STERILE PVC TUBING ON REEL			
CODE	DIAMETER TUBE	LENGTH	PACKAGING
H93005N	8 X 12 mm	50 m	2 reels
H93006N	10 X 14 mm	40 m	2 reels

