


URINE BAGS

ECONOMIC RANGE



- Good quality at attractive price.
- Different options ; with handle, Pasteur drip chamber.
- Latex free.



BioService
TUNISIE
Quality is not a coincidence ...



ECODRAIN:

- The « Ecodrain » is the optimal choice for an Economic Drainage. It offers a transparent good quality of urine bag at very attractive price.
- Available in 2 options, with and without handle.
- BioService Tunisie is able to supply it in different tube length according to user's requirements.

ECODRAIN			
REFERENCE	MODEL	TUBE LENGTH	PACKAGING
H90240E	ECODRAIN 2000-90	90 cm	100 pcs
H90241E	ECODRAIN 2000-100	100 cm	100 pcs
H90242E	ECODRAIN 2000-110	110 cm	100 pcs
H90250E	ECODRAIN 2000-130	130 cm	100 pcs

Your safety is our priority

BIO SAFETY:

- Since the safety of our patient's is quite important, BioService Tunisie thought about all possible requirements to grant you new models of urine bags with different specifications.
- « BioSafety » contains the necessary components for a safety urine collection, with Pasteur Drip Chamber, Non return valve, needle-free connector for easy sample collection without needle and drainage T-tap with secure tap support.
- Available in 3 models :

BIO SAFETY			
REFERENCE	MODEL	TUBE LENGTH	PACKAGING
H90335E	BioSafety 2000-120	120 cm	50 pcs
H90465E	BioSafety 2000-120 with handle	120 cm	40 pcs
H90565E	BioSafety 2000-120 with Pasteur drip chamber	120 cm	40 pcs

The bags are available in standard conformation in individual packaging, and ETO sterilized